

SPACE BRACE

Merkblatt



Hauptgeschäft Weitere Standorte

Effingerstrasse 37 CH-3008 Bern Fon 031 388 89 89 Fax 031 382 34 81	Bälliz 55 CH-3600 Thun Fon 033 222 30 55 Fax 033 222 30 26	Hammerstrasse 12 CH-4058 Basel Fon 061 691 62 70 Fax 061 693 15 52	Aaraustrasse 42 CH-5734 Reinach AG Fon 062 772 36 60 Fax 062 772 36 61	Forchstrasse 313 CH-8008 Zürich Fon 044 422 46 22 Fax 044 422 46 30	Oberer Graben 28 CH-8400 Winterthur Fon 052 209 01 90 Fax 052 209 01 91
--	---	---	---	--	--

www.ortho-team.ch
info@ortho-team.ch

SPACE BRACE



1 Klettverschluss oben schliessen

2 Kreuzzuge schliessen

3 Klettverschluss Sprunggelenk schliessen

4 Ballen Velcro schliessen

Bitte beachten! Fuss muss beim Anziehen im rechten Winkel stehen.